



XENOMETRIX

Ames II


Ames MPF™ 98/100

Ames MPF™ 1535

Ames MPF™ 1537

Ames MPF™ E.coli


Ames MPF™ PENTA I



- 
- The Ames MPF and Ames II Fluctuation Assays come as complete ready to use kits with strains, liquid culture media, positive control and S9
  - The Ames MPF and Ames II Fluctuation Assays can be used to detect genotoxic activity in chemicals, cosmetics, pharmaceuticals, foodstuffs, water, air, soils or sediments.
  - Quality controlled *S. typhimurium* and *E.coli* strains (genotyped and phenotyped)
  - 3 times less compound required
  - 5 times less operator intervention
  - 3 fold less contaminated waste
  - Less consumption of rat liver S9
  - 10 times lower cost for analysis
  - OECD Guideline 471 and FDA compliant
  - Colorimetric determination to reduce reporting errors

# Ames MPF, Ames II

## Test procedure



### Please send me more information:

- Instruction Manual:     Ames MPF 98                       Ames MPF 100                       Ames MPF 1535
- Ames MPF 1537                       Ames MPF E.coli                       Ames MPF 98/100
- Ames MPF PENTA I (TA98/100/1535/1537/E.coli)                       Ames II
- Price list
- References
- Please contact me

Name: \_\_\_\_\_ First Name: \_\_\_\_\_

Affiliation: \_\_\_\_\_

Street: \_\_\_\_\_

Zip/City/Country: \_\_\_\_\_

Email/Tel.: \_\_\_\_\_

Xenometrix AG  
 Gewerbestrasse 25  
 CH-4123 Allschwil  
 Switzerland

Tel ++41-61-482 14 34  
 Fax ++41-61-482 20 72  
 info@xenometrix.ch  
 www.xenometrix.ch

Your local distributor: